

















# AESCULAP Aicon® Instrument Organisation System

## Assignment of instrument categories to instrument holders

The following overview provides comprehensive understanding of the assignment of instrument categories to the different instrument holders of the AESCULAP Aicon® Instrument Organization System. The assignment is only a suggestion and not binding. Please always follow the IFU of both, the instruments and the instrument holders. If you would like a personal consultation, kindly contact us.

																
	JJ910	JJ911	JJ912	JJ913	JJ914	JJ915	JJ916	JJ917	JJ918	JJ919	JJ920	JJ921	JJ922	JJ923	JJ924	JJ925
Bone levers			✓	✓												
Cannulas			✓	✓	✓	✓										
Chisels							✓	✓								
Clamps			✓	✓	✓	✓										
Curettes			✓	✓				✓	✓	✓	✓	✓	✓		✓	
Dilators			✓	✓					✓	✓	✓	✓	✓		✓	
Dissectors			✓	✓												
Electrodes / Electrode handles									✓	✓	✓	✓	✓		✓	
Elevators							✓	✓								
Files							✓	✓								
Forceps (Tweezer)			✓	✓												
Hooks			✓	✓												
Instruments with round shafts of various diameters									✓	✓	✓	✓	✓		✓	
Instruments with sliding or tubular shaft			✓	✓	✓	✓										
Mallets								✓	✓	✓	✓	✓	✓		✓	
Needle Holders			✓	✓												
Osteotomes							✓	✓								
Probes					✓	✓										
Punches			✓	✓												

No special instrument group.  
Dividers are intended to separate e. g. heavy instruments such as mallets from delicate instruments or to categorize different instrument groups.

# AESCULAP Aicon® Instrument Organisation System

Assignment of instrument categories to instrument holders



Raspatory						✓	✓									
Redon guide needles			✓	✓	✓	✓										
Rongeurs			✓	✓												
Scalpel handles			✓	✓	✓	✓										
Scissors			✓	✓												
Scopes	✓	✓						✓	✓	✓	✓	✓	✓	✓		
Strippers					✓	✓										
Suction Instruments			✓	✓	✓	✓										

No special instrument group.  
Dividers are intended to separate e. g. heavy instruments such as mallets from delicate instruments or to categorize different instrument groups.


## Examples



JJ912



JJ912 / JJ913



JJ914 / JJ915




JJ917




JJ918 / JJ919 / JJ920 / JJ921 / JJ922 / JJ923




## Combination possibilities for fixation clamps & instrument holders




AESCULAP Aicon® fixation clamp small (JJ901R) can be combined with one small instrument holder (e. g. JJ918)



AESCULAP Aicon® fixation clamp medium (JJ902R) can be combined with one medium instrument holder (e. g. JJ917)



AESCULAP Aicon® fixation clamp large (JJ903R) can be combined with four small instrument holders (e. g. JJ912, JJ920, JJ919, JJ916) or with two medium instrument holders (e. g. JJ913, JJ923) or with one medium instrument holder (e. g. JJ923) and two small instrument holders (e. g. JJ920, JJ924)



For mounting instructions please scan the QR-code or visit the website: <https://bbraun.info/aos>